



BIOSCANNER™ TOOLS

ADAPT OVERCOME PROSPER



SCANNING TABLE

Solid wood table with acrylic centerpiece facilitate of scans that require supine positioning. Clear straps assist with the proper positioning of the patient. May be used for TLSO, LSO or CTLSO. Will work with BioScanner which has scan-through-glass technology to correct for refraction.

SCANNING PANTS

Black and White knit pants designed specifically for scanning. Black area allows scans to be taken effortlessly without interference from the contralateral leg. May be used for TF, KAFO and HKAFO.

CRANIAL SCANNING HOOD

Flattens hair and allows user to see prominence of ears and neck to create an accurate fitting cranial helmet.

