



Kinetec* *Maestra* Hand & Wrist CPM Machine

The
CPM
concept



Rehabilitation

Post-Operatively



- **The CPM Solution**
handling every single hand and wrist pathology
- **A useful tool for protocols**
dedicated to surgeons, physiotherapists
and occupational therapists
- **Rehabilitation/Post-operatively**
- **Each hand and wrist movement possible**
- **Analytic / Ambulatory**
- **Kinetec Quality**



The protocols' tool which helps in finding the CPM solution

CPM SOLUTION


POST-OPERATIVE

Maestra Portable

3 Movements
Long Fingers



1• Composite fist
(MCP-PIP-DIP)
-15° to 270°



2• Intrinsic -
0° to 180°



3• Intrinsic +
0° to 90°



Ordering Information

Complete apparatus with pads 46210055

INDICATIONS

- Flexor and extensor tenolysis.
- Aponeurectomies for Dupuytren's disease.
- Open reduction and internal fixation of intra-articular, diaphyseal, metaphyseal and epiphyseal fractures of the phalanges.
- MCP arthrolysis.
- Prosthetic replacement of the MCP, PIP and DIP joints.
- Rheumatoid/neurological or after-burn stiffness.
- Capsulotomies, arthrolysis and tenolysis for post-traumatic stiffness of MCP, PIP and DIP joints.
- Synovectomies.

REHABILITATION

Maestra

9 Movements
Hand & Wrist



1• Composite fist
(MCP-PIP-DIP)
-30° to 225°



- 2• Ulnar deviation
-30° to 60°
- 3• Extension / Flexion of the wrist
-30° to 60°
- 4• Extension / Flexion of the wrist
with MCP
-50° to 140°
- 5• Pronation / supination
-90° to 90°



6• Thumb Opposition
0° to 180°



7• Intrinsic +
0° to 90°

8• Intrinsic -
0° to 180°

9• Selective mobilisation
of the DIP joint
0° to 70°



Ordering Information

Complete apparatus with pads 46210050

KINETEC

AbilityOne

AbilityOne Kinetec S.A.
Rue Albert Deville - Tournes
F-08014 Charleville-Mézières Cedex
Phone : +33 (0) 3.24.52.91.21
Fax : +33 (0) 3.24.52.90.34

Your local distributor is :



Sold & Serviced In India by :

BR Biomedicals

C-2/2536, Vasant Kunj, New Delhi-110 070

Tel.: 011-2689 7274 Fax: 011 - 2689 8757

E-mail : sales@brbiomedicals.com

Website : www.brbiomedicals.com

* Kinetec is a trademark of AbilityOne.

AbilityOne Kinetec reserves the right to make modifications not described in this brochure.

4696710184 - 11/2004